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## Instructions for Use

### Crown Temp (AUTOMIX 10:1)

#### Automix Temporary Crown & Bridge Material

Crown Temp (AUTOMIX 10:1) is an **automatically mixed** two component material **in cartridges** based on multifunctional methacrylic esters. The temporary material is a very easy and comfortable material for creating **temporary crowns or bridges**. Because of its flexibility and high hardness, the material is especially suitable for longer bridge spans. Crown Temp (AUTOMIX 10:1) shows fluorescence in UV-light.

Crown Temp (AUTOMIX 10:1) is free of methyl methacrylate. Its temperature derived while curing is lower than 40 °C. As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows increasing transverse strength, good abrasion resistance and low polymerization shrinkage. It fits perfect. Showing good polishability, good color stability and fluorescence it looks perfect.

And if Crown Temp (AUTOMIX 10:1) breaks - which may happen with acrylics - it is easy to repair. Cured material can be repaired with any natural or light cure composite building up again a chemical link. In nearly every property Crown Temp (AUTOMIX 10:1) is superior to the common materials based on monomer/polymer systems. It allows thin edges and production of longer bridge spans.

#### Important advantages of Crown Temp (AUTOMIX 10:1):

- no dosage and mixing by hand, no transfer into an application syringe, exact dosage,
- direct bubble-free application
- no extra manipulation time necessary
- excellent fit (small polymerisation shrinkage)
- excellent aesthetics (polishability, translucency, colour stability, fluorescence)
- very good durability

#### Application:

##### 1. Impression taking

Before preparing of stumps or before a planned extraction of a tooth a situation impression is made by addition curing silicones (long storage stability!) or by alginate. Interdental areas are carved out, in molar areas with teeth absent it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

##### 2. Preparing of the cartridge

The turn cap of the safety cartridge is taken off and thrown away (do not use it again!). It is substituted by the special Crown Temp (AUTOMIX 10:1) mixing cannula by attaching it with a sideways turning movement, locking the cannula in place. The application gun is loaded with the prepared cartridge and is ready for application.

###### Note:

The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.

Leave the used mixing cannula on the cartridge. It serves as a cap.

##### 3. Application

Crown Temp (AUTOMIX 10:1) is automatically mixed when dispensed with slight and even pressure directly into the situation impression made before. Filling should occur from bottom upward to prevent voids.

##### 4. Forming of the temporary crowns or bridges

- a. Load the situation impression with Crown Temp .
- b. Seat the impression onto the prepared areas of the teeth. Step a and b must be done within 60 seconds (working time).
- c. after 1-2 minutes (setting time in mouth) the material shows a hardened but still elastic condition and can be removed from the teeth (alternatively: 3-4 minutes setting time on the model)

###### Note:

The temporary crown and bridge substitute can only be removed without destruction during the elastic state

## 5. Post curing and finishing

After removal of the temporary crown or bridge from the impression (or from the tooth preparation) excess material and proximal undercuts are removed. Preferably, the temporary is then post cured in warm (45°C - 55°C) water (e.g. in a hot cure polymerization device). If this is not possible or not desired, after **6 minutes from beginning of application** the temporary acrylic can be worked out with rotary instruments and can be polished to high gloss.

###### Note:

The oxygen inhibited smear layer on the surface caused by the acrylic system can not be fully avoided and should be removed before working out. It can easily be removed by alcohol or other suitable solvents.

## 6. Cementing of the temporary

Crown Temp (AUTOMIX 10:1) temporary crowns or bridges preferably should be cemented with a **eugenol-free temporary cement** (e.g. Crystal Temp NE). In case of using eugenol-containing cements it should be considered that later used acrylate based cements could be hindered in hardening (inhibition of curing reaction by traces of remained eugenol).

## 7. Repairs

Crown Temp (AUTOMIX 10:1) temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

### a) Fracture shortly after production

Both ends of the fracture are refixed with freshly extruded Crown Temp (AUTOMIX 10:1).

### b) Fracture of longer existing temporary

The areas of fracture are cleaned and roughened and provided with some mechanical retentions. Thus prepared, the fractured areas can be joint with freshly extruded Crown Temp (AUTOMIX 10:1) To enhance the curing place for some minutes in warm water.

## 8. Special hints

Non hardened Crown Temp (AUTOMIX 10:1) can easily be removed with alcohol or other suitable solvents.

## 9. Precautions

Crown Temp (AUTOMIX 10:1) is free of methyl methacrylate but contains other methacrylates.

With susceptible patients, sensitization to Crown Temp (AUTOMIX 10:1) cannot be excluded. In these cases Crown Temp (AUTOMIX 10:1) should not be used any more, if allergic reactions are observed.

Do not use for patients with allergic reactions against acrylates.

## 10. Storage

Do not store above 20 °C!

Crown Temp (AUTOMIX 10:1) shall not be used after expiry date!

**Contents of cartridge: 50 ml**